MIURA Technical Report

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Rapid Detection Method of Viable *Legionella* Species by the Combination of Continuous-Flow Dielectrophoresis System and Quantitative PCR

Yuko MIYAUCHI¹, Honami IKEUCHI¹, Masaki TAKAI¹, Ayato KAWASHIMA²,
Satoshi MORITA³, Yusei SATO³, and Takaharu ENJOJI³

¹ MIURA CO., LTD., 864-1 Hojotsuji, Matsuyama, Ehime 799-2430, Japan

² Graduate School of Agriculture, Ehime University, 3-5-7 Tarumi, Matsuyama, Ehime 790-8566, Japan

³ AFI Corporation, 2nd Flr., Medical Innovation Center, Kyoto University,
53 Shogoin Kawahara-cho, Sakyo-ku, Kyoto 606-8397, Japan

Summary

In this study, we investigated the method of separating and quantifying only viable *Legionella* by combining the Continuous-Flow Dielectrophoresis System (CFDS) and quantitative PCR (qPCR). The CFDS is equipped with a particulate filtration system that combines the dielectrophoresis (DEP) principle and microfluidic control technology. The DEP force has a frequency dependence, and when the separation state of viable and dead

Legionella by CFDS was confirmed with a CCD camera, the capture of viable *Legionella* was observed significantly in a specific frequency range (3.0 MHz to 6.0 MHz). So, the trapped rate of viable *Legionella* at an applied voltage of 22.0 Vpp and frequency of 6.0 MHz was determined. Furthermore, the number of copies of 16S rRNA gene per 1 CFU was calculated as a conversion factor based on the number of colonies by the plate culture method and the number of copies by qPCR from *Legionella* recovered samples captured by CFDS (22.0 Vpp, 6.0 MHz). The conversion formula was then established to determine the number of *Legionella* by this method. Samples prepared at nine different bacteria concentrations were then analyzed for the number of *Legionella* by this method and compared with the plate culture method. As a result, $R^2 = 0.9353$ was obtained, indicating a high correlation between this method and the plate culture method.

Key words: Dielectrophoresis (誘電泳動) / *Legionella* Species (レジオネラ属菌) / Rapid detection method (迅速検出法) / Quantitative PCR (qPCR)

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
Fig. 1. Setting up for analysis using Continuous-Flow Dielectrophoresis System (CFDS).

A: Hardware of CFDS “ELESTA PixeeMo”, B: Component parts – a) Syringe pump unit, b) Digital Microscope, c) “ELESTA Chip” for separating viable and dead bacteria, d) Syringe for pumping or reverse pumping liquid, e) PFA tube, f) Reverse liquid flow unit, C: Configuration of “ELESTA Chip” – g) Inlet, h) Electrodes for DEP, i) Outlet, D: Magnified image of electrodes by CCD camera – j) Electrode, k) Slit between electrodes, l) Bacteria trapped by CFDS.

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

Fig. 2. Images of viable or dead *Legionella* trapped by CFDS at each frequency. Quantity level of trapped bacteria: (++) high, (+) middle, (-) low.

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Table 1. Trapped rate of viable *Legionella* by CFDS at 3.0 MHz or 6.0 MHz (N=3).

Frequency (MHz)	3.0	6.0
Trapped rate of viable <i>Legionella</i> (%)	42 $\pm$ 2.5	31 $\pm$ 2.3

Table 2. Residual rate of dead *Legionella* by CFDS at 3.0 MHz or 6.0 MHz (N=3).

Frequency (MHz)	3.0	6.0
Residual rate of dead <i>Legionella</i> (%)	0.096 $\pm$ 0.055	0.015 $\pm$ 0.0034

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Table 3. Detection results of *Legionella* by the plate culture or qPCR in each 100 μL of recovery sample by CFDS and the conversion factor calculated from 9 combinations.

Initial bacteria concentration (CFU/mL)	Plate culture method (CFU/100 μL)		qPCR (copies/100 μL)		Conversion factor (copies/CFU)					
10^2	Sample1 (S1)	11	Sample4 (S4)	75.1	S4/S1	6.8	S5/S1	5.9	S6/S1	6.0
	Sample2 (S2)	12	Sample5 (S5)	64.9	S4/S2	6.3	S5/S2	5.4	S6/S2	5.5
	Sample3 (S3)	14	Sample6 (S6)	65.9	S4/S3	5.4	S5/S3	4.6	S6/S3	4.7
					Avg.=5.6 SD=0.71 CV(%)=13					
10^3	Sample7 (S7)	101	Sample10 (S10)	514	S10/S7	5.1	S11/S7	5.8	S12/S7	6.3
	Sample8 (S8)	93	Sample11 (S11)	585	S10/S8	5.5	S11/S8	6.3	S12/S8	6.8
	Sample9 (S9)	102	Sample12 (S12)	635	S10/S9	5.0	S11/S9	5.7	S12/S9	6.2
					Avg.=5.9 SD=0.60 CV(%)=10					

$$\text{元試料菌数 [CFU/mL]} = \frac{\text{定量コピー数 [copies/回収液]}}{6} \div 0.3 \quad (\text{b})$$

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
Fig. 3. Comparison of the detection results of *Legionella* by the plate culture and CFDS + qPCR. The result of CFDS + qPCR is applied to the conversion formula.

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